

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

 β -Glucuronidase from Helix pomatia - Type H-3, aqueous solution, $\geq 90,000$ units/mL

Product Number:

G8885

MDL:

MFCD00130629

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

unit/ml β -Glucuronidase ≥ 90000

One unit will liberate 1.0 Microgram of Phenolphthalein from Phenolphthalein glucuronide per hour at pH 5.0 at 37 Deg C.

Sulfatase units/ml

 ≤ 1000

One unit will hydrolyze 1.0 Micromole of p-Nitrocatechol sulfate per hour at pH 5.0 at 37 deg C

Specification: PRD.1.ZQ5.10000034660

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.