


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1-Aminoindan - 98%

Product Number: A59506
CAS Number: 34698-41-4
MDL: MFCD00003799
Formula: C₉H₁₁N
Formula Weight: 133.19 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
------	---------------

Appearance (Color)	Faint Green-Yellow to Green-Yellow and Faint Yellow to Dark Yellow
Appearance (Form)	Liquid
Proton NMR spectrum	Conforms to Structure
Purity (GC)	≥ 97.5 %

Specification: PRD.1.ZQ5.10000045833

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.