

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium chloride - Vetec™ reagent grade, 99%

Product Number:

V900068

CAS Number:

7447-40-7

Formula:

KCl

Formula Weight:

74.55 g/mol

KCl

TEST	Specification
------	---------------

Appearance (Colour)

White

Appearance (Form)

Powder

Solubility (Colour)

Colorless

Solubility (Turbidity)

Clear

20 g + 150 mL, H₂OTitration by AgNO₃

≥ 99 %

Specification: PRE.0.ZQ5.10000018410

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.