

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2,2-Dimethyl-3,6,9,12-tetraoxa-2-silatridecane - ≥98%, anhydrous

Product Number:

900871

CAS Number:

864079-62-9

Formula:

C₁₀H₂₄O₄Si

Formula Weight:

236.38 g/mol

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared Spectrum	Conforms to Structure
¹ H NMR Spectrum	Conforms to Structure
Purity (GC)	≥ 98.0 %
Water (by Karl Fischer)	≤ 500 ppm
Boiling point by DSC	
Specification: Report Only	

Specification: PRD.0.ZQ5.10000088327

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.