


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Triton™ X-45

Product Number: **X45**
CAS Number: 9002-93-1
MDL: MFCD00128254
Formula: $(C_2H_4O)_n C_{14}H_{22}O$



TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Viscous Liquid
Brookfield Viscosity	260 - 370 cps
c = Neat At 25 Degrees Celsius	
Cloud Point	39.4 - 44.8 °C
Registered Trademark	Confirmed
Triton is a trademark of The Dow Chemical Company or an affiliated company of Dow	

Specification: PRD.1.ZQ5.10000030259

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.