

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Undecanone - ≥98%, FCC, FG

Product Number: W309303
CAS Number: 112-12-9
MDL: MFCD00009583
Formula: C₁₁H₂₂O
Formula Weight: 170.29 g/mol



TEST	Specification
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.428 - 1.432
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.0 %
Specific Gravity At 25°C	0.822 - 0.826
Acid Value	≤ 5.0 ml
Solubility (Turbidity) 1ml/1ml 95% Ethanol	Clear
Solubility (Color)	Colorless to Faint Yellow
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period 5 Years	-----

Specification: PRD.3.ZQ5.10000005842

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.