

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platinum on carbon - extent of labeling: 10 wt. % loading, matrix activated carbon support

Product Number: 205958

MDL: MFCD00011179

Formula: Pt

Formula Weight: 195.08 g/mol

Pt

TEST	Specification
------	---------------

Appearance (Color)

Black

Appearance (Form)

Powder

Purity

% Platinum

9.8 - 10.2 %

Loss on Drying

≤ 3.0 %

Specification: PRD.0.ZQ5.10000076628

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.