

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Puromycin dihydrochloride - Ready Made Solution, from Streptomyces alboniger, 10 mg/mL in H₂O, 0.2 µm filtered

Product Number:


P9620

Formula Weight:

544.44 g/mol

Storage Temperature:

-20 °C



TEST	Specification
------	---------------

Appearance (Form)

Clear Solution

Conforms

Appearance (Colour)

Colorless to Faint Yellow

Concentration

10.0 - 10.2 mg/ml

Purity (HPLC)

≥ 98.00 %

Cell Culture (Suitability)

Pass

Recommended Retest Period

2 Years

Specification: PRE.0.ZQ5.10000010047

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.