

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
SIGMAFAST™ Protease Inhibitor Tablets - For General Use

Product Number: **S8820**

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	White to Off White
Appearance (Form)	Round Tablet
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear to Hazy
	One tablet dissolved in 100 mL of H ₂ O with gentle agitation gives a clear to hazy colorless solution.
Suitability	Pass
	A 10X concentrate solution of Inhibitor will inhibit a 1.62 mg/ml suspension of P1625 > or = 95%

Specification: PRD.2.ZQ5.10000027344

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.