


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Chitosan - low molecular weight

Product Number: **448869**  
CAS Number: 9012-76-4  
MDL: MFCD00161512



TEST	Specification
Appearance (Color)	Faint Beige to Beige
Appearance (Form)	Conforms to Requirements
Powder and/or Chips	
Assay	≥ 75 %
Deacetylation	
Viscosity	20 - 300 cps
c = 1% in 1% Acetic Acid	

Specification: PRD.2.ZQ5.10000020770

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.