


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Thebaine - powder

Product Number: **T2019**
CAS Number: 115-37-7
MDL: MFCD00027409
Formula: C₁₉H₂₁NO₃
Formula Weight: 311.37 g/mol



TEST	Specification
Solubility (Turbidity)	Clear
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
¹ H NMR Spectrum	Conforms to Structure
Titration	≥ 98 %
Purity (HPLC)	≥ 98 %

Specification: PRD.2.ZQ5.10000048657

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.