


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Oil Red EGN

Product Number: **234117**
CAS Number: 4477-79-6
MDL: MFCD00009688
Formula: C₂₅H₂₂N₄O
Formula Weight: 394.47 g/mol



TEST	Specification
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Wavelength c = 0.02g/L; Chloroform	517 - 525 nm
Extinction Coefficient	≥ 25000
Wavelength	352 - 360 nm
Extinction Coefficient	≥ 14000
Carbon	≥ 64.7 %
Nitrogen	≥ 12.1 %

Specification: PRD.0.ZQ5.10000035396

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.