


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Pentylentetrazole

Product Number: P6500
CAS Number: 54-95-5
MDL: MFCD00005939
Formula: C₆H₁₀N₄
Formula Weight: 138.17 g/mol
Storage Temperature: -20 °C



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
50 mg/mL, H ₂ O	
Carbon	51.4 - 52.9 %
Nitrogen	40.0 - 41.1 %
Purity (GC)	≥ 99 %

Specification: PRD.2.ZQ5.10000008004

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.