

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

 $\beta$ -Glucuronidase from abalone - aqueous solution,  $\beta$ -glucuronidase  $\geq 100,000$  units/mL

Product Number:

SRE0037

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

unit/ml  $\beta$ -Glucuronidase $\geq 100000$ 

Unit Definition: One unit will liberate 1.0 microgram of phenolphthalein from phenolphthalein glucuronide per hour at pH 3.8 at 37 deg C. (30 minute assay)

Sulfatase units/ml

 $\leq 8000$ 

Unit Definition: One unit will hydrolyze 1.0 micromole of p-nitrocatechol sulfate per hour at pH 5.0 at 37 deg C (30 minute assay).

Specification: PRD.1.ZQ5.10000056533

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.