

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2-Chloro-4-fluoro-5-sulfamoylbenzoic acid - 96%

Product Number:

479845

CAS Number:

4793-24-2

MDL:


MFCD01321357

Formula:

C₇H₅ClFNO₄S

Formula Weight:

253.64 g/mol

**TEST****Specification**

Appearance (Color)

White

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Titration by AgNO₃

95.5 - 104.5 %

% Purity Based On Chloride Content By Titration With AGNO₃ After O₂

Combustion

Purity (HPLC)

≥ 95.5 %

Specification: PRD.0.ZQ5.10000076505

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.