

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2-Ethoxycarbonylphenylboronic acid - ≥95%

Product Number:

683833

CAS Number:

380430-53-5

MDL:


MFCD02179453

Formula:

C₉H₁₁BO₄

Formula Weight:

193.99 g/mol

**TEST****Specification**

Appearance (Color)

White to Tan

Appearance (Form)

Powder or Crystals

Proton NMR Spectrum

Conforms to Structure

Purity (GC)

≥ 95.0 %

Note

Note

May Contain Varying Amounts of Anhydride

Specification: PRD.1.ZQ5.10000036528

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.