

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
2,2'-Azobis(2-methylpropionitrile) - recrystallized from methanol, 99%

Product Number: **755745**

CAS Number: 78-67-1



Formula Weight: 164.21 g/mol

Storage Temperature: -20 °C

### TEST

### Specification

Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Proton NMR spectra	Conforms to Structure
Carbon	57.6 - 59.4 %
Nitrogen	33.6 - 34.6 %

Specification: PRD.2.ZQ5.10000022032

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.