


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Methyl 3-oxotetrahydrothiophene-2-carboxylate - 96%

Product Number: 650633
CAS Number: 2689-69-2
MDL: MFCD00156156
Formula: C₆H₈O₃S
Formula Weight: 160.19 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Carbon	43.0 - 47.0 %
Sulfur	19.1 - 20.9 %
Purity (GC)	≥ 95.5 %

Specification: PRD.0.ZQ5.10000029491

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.