


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Cellulose - colloidal, microcrystalline

Product Number: **435244**
CAS Number: 9004-34-6
MDL: MFCD00081512



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Sodium Carboxymethyl Cellulose	10.0 - 20.0 %
Loss on Drying	≤ 8.0 %
Size: + 60 Mesh	≤ 0.10 % wt
Size: + 325 Mesh	≤ 70.0 % wt
Heavy Metal	≤ 0.001 %

Remarks:

Reviewing Chemist :

Specification: PRD.0.ZQ5.10000101278

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.