


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Chiralyst P493 - Umicore, 98%

Product Number: **683132**
CAS Number: 12092-47-6
MDL: MFCD00012415
Formula: C₁₆H₂₄Cl₂Rh₂
Formula Weight: 493.08 g/mol



TEST

Specification

Appearance (Color)
Orange, Yellow-Orange or Orange-Brown
Appearance (Form)
Infrared spectrum
Rhodium (Rh)
Note
Product of Umicore

Conforms to Requirements

Powder or Crystals
Conforms to Structure
41.0 - 42.0 %
Note

Specification: PRD.1.ZQ5.10000012339

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

