

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
SULFUR DIOXIDE SOLUTION, 0.5 M IN THF

Product Number: 901592

TEST	Specification
Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
ICP Major Analysis Confirms Sulfur Component	Confirmed
¹ H NMR Spectrum THF	Conforms to Structure
Titration with Na ₂ S ₂ O ₃	0.45 - 0.55 M

Specification: PRD.0.ZQ5.10000088297

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.