

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Gelatin from porcine skin - gel strength ~175 g Bloom, Type A

Product Number:

**G2625**

CAS Number:

9000-70-8

MDL:

MFCD00081638

**TEST****Specification**

Appearance (Color)

Faint Yellow to Light Yellow to Light Beige

Appearance (Form)

Powder

Solubility (Color)

Faint Yellow to Yellow

Solubility (Turbidity)

Slightly Hazy to Very Hazy

50 mg/mL, H<sub>2</sub>O

Protein by Biuret

65 - 85 %

Gel Strength

160 - 190

(Bloom Gelometer)

(approx 175)

Specification: PRD.2.ZQ5.10000026360

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.