

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

4,4'-(Oxydi-2,1-ethanediyl)bismorpholine - 97%

Product Number:

424528

CAS Number:

6425-39-4

MDL:


MFCD00072740

Formula:

C₁₂H₂₄N₂O₃

Formula Weight:

244.33 g/mol

**TEST****Specification**

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (TLC)

≥ 97 %

Solubility (Turbidity)

Clear

50mg/ml, Methanol

Solubility (Color)

Colorless to Yellow

Specification: PRD.0.ZQ5.10000050984

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.