

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Phenylphosphonic dichloride - technical grade, 90%

Product Number:

389560

CAS Number:

824-72-6

MDL:


MFCD00002070

Formula:

C₆H₅Cl₂OP

Formula Weight:

194.98 g/mol

**TEST****Specification**

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 88.0 %

Specification: PRD.0.ZQ5.10000028984

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.