

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Activated charcoal - DARCO®, 4-12 mesh particle size, granular

Product Number: 242233
CAS Number: 7440-44-0
MDL: MFCD00133992
Formula: C
Formula Weight: 12.01 g/mol

C

TEST	Specification
------	---------------

Appearance (Color)	Black
Appearance (Form) Small Chunks	Conforms
Particle Size Distribution >4 Mesh	≤ 5 %
Particle Size Distribution <12 Mesh	≤ 5 %
Moisture Content	≤ 12 %
Dust	≤ 0.40 %
Registered Trademark DARCO is a registered trademark of Norit Americas, Inc.	-----

Remarks:

Reviewing Chemist :

Specification: PRD.0.ZQ5.1000028161

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.