


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Succinic acid disodium salt - 99%

Product Number: 224731
CAS Number: 150-90-3
MDL: MFCD00002790
Formula: C₄H₄Na₂O₄
Formula Weight: 162,05 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Titration with HClO ₄	98.5 - 101.5 %

Specification: PRD.1.ZQ5.10000019956

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.