

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Triethylenetetramine - technical grade, 60%

Product Number: 132098
CAS Number: 112-24-3
MDL: MFCD00008169
Formula: C₆H₁₈N₄
Formula Weight: 146.23 g/mol



TEST	Specification
Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Liquid
GCMS	Conforms to Structure
Purity (GC)	≥ 59.5 %
Color Test	≤ 150 APHA

Specification: PRD.2.ZQ5.10000008673

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.