

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Aluminum - $\leq 30 \mu\text{m}$, spherical powder, $\geq 99\%$ **Product Number:****266515**

Al

CAS Number:

7429-90-5

MDL:

MFCD00134029

Formula:

Al

Formula Weight:

26.98 g/mol

TEST**Specification**

Appearance (Color)

Conforms to Requirements

Grey

Appearance (Form)

Powder

Complexiometric EDTA % Al

99.0 - 101.0 %

X-Ray Diffraction

Conforms to Structure

Size

 ≤ 30 micron

Average Particle Size

Specification: PRD.0.ZQ5.10000026465

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.