

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Polyethylene - Ultra-high molecular weight, surface-modified, powder, 125 µm avg. part. size

Product Number:

434264


CAS Number:

9002-88-4

MDL:

MFCD00084423

Formula:

C₂H₄

TEST	Specification
------	---------------

Appearance (Color)

White

Appearance (Form)

Conforms to Requirements

Fine Powder

Infrared spectrum

Conforms to Structure

Size

105 - 145 micron

Specification: PRD.1.ZQ5.10000060756

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.