


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ammonium acetate - 99.999% trace metals basis

Product Number: 372331
CAS Number: 631-61-8
MDL: MFCD00013066
Formula: C₂H₇NO₂
Formula Weight: 77.08 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Crystals
Infrared spectrum	Conforms to Structure
Titration by HCL	≥ 97.0 %
Trace Metal Analysis	≤ 20.0 ppm
Purity	Conforms
99.999% Based on Trace Metals Analysis	

Specification: PRD.0.ZQ5.10000026144

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.