

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Potassium iodide - ReagentPlus®, 99%

**Product Number:** 207969  
**CAS Number:** 7681-11-0  
**MDL:** MFCD00011405  
**Formula:** IK  
**Formula Weight:** 166.00 g/mol

**KI**

<b>TEST</b>	<b>Specification</b>
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
X-Ray Diffraction	Conforms to Structure
Titration with AgNO <sub>3</sub> % Iodine	75.3 - 77.6 %
Recommended Retest Period 42 Months	-----

Specification: PRD.3.ZQ5.10000016864

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.