


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium selenite – BioReagent, suitable for cell culture, ≥98%

Product Number: S5261
CAS Number: 10102-18-8
MDL: MFCD00003489
Formula: Na₂O₃Se
Formula Weight: 172.94 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear to Very Slightly Hazy
50 mg/mL, H ₂ O	
Identity - Selenium (Se) (ICP)	Confirmed
Purity by Titration	≥ 98 %
Redox Titration	
Cell Culture Test	Pass
Recommended Retest Period	-----
4 Years	

Specification: PRD.0.ZQ5.10000027656

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.