

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Osteoprotegerin antibody produced in goat – affinity isolated antibody, lyophilized powder

Product Number: **O1139**

Storage Temperature: -20 °C

TEST	Specification
Formulation Lyophilized from 0.2 um-filtered solution in phosphate buffered saline containing carbohydrates	Conforms
Endotoxin Assay < 0.1 EU/ug antibody	Conforms
Bioactivity ND50: 0.5 - 2 ug/mL	Conforms
Working Dilution by Indirect Elisa 0.5-1 ug/mL detects approximately 0.5 ng/well of recombinant mouse OPG	Conforms
Working Dilution by Indirect Immunoblotting 0.1-0.2 ug/mL detects 2 ng/lane under non-reducing and 5 ng/lane under reducing conditions	Conforms

Specification: PRD.0.ZQ5.10000070485

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.