

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Polypropylene, methyl-PSS nanoreinforced - contains ~10 wt. % methyl-PSS as reinforcer

Product Number:

565628

MDL:

MFCD00084447



TEST

Specification

Appearance (Color)

White to Off-White

Appearance (Form)

Beads

Infrared spectrum

Conforms to Structure

Size

Specification: Report Result. Units: mm

Specification: PRD.0.ZQ5.10000024730

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.