

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Vitamin K1 - viscous liquid


Product Number: **V3501**

CAS Number: 84-80-0

Formula: C₃₁H₄₆O₂

Formula Weight: 450.70 g/mol

Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)	Yellow to Orange
Appearance (Form)	Viscous Liquid
Solubility (Color)	Yellow to Dark Yellow
Solubility (Turbidity)	Clear
100 ul/ml, CHCl ₃	
Infrared spectrum	Conforms to Structure
Carbon	81.0 - 84.3 %
Recommended Retest Period	-----
2 years	

Specification: PRD.0.ZQ5.10000018220

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.