

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Calcium sulfate -  $\geq 99.99\%$  trace metals basis**Product Number:****255696** $\text{CaSO}_4$ 

CAS Number:

7778-18-9

MDL:

MFCD00010912

Formula:

CaO4S

Formula Weight:

136.14 g/mol

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Chunks
Complexiometric EDTA % Calcium	28.5 - 30.3 %
ICP Major Analysis Confirms Calcium Component	Confirmed
Trace Metal Analysis	$\leq 100.0$ ppm
Purity $\geq 99.99\%$ Based On Trace Metals Analysis	Meets Requirements

Specification: PRD.0.ZQ5.10000025606

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.