

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Lithium sulfate -  $\geq 99.99\%$  trace metals basis**Product Number:****203653**

CAS Number:

10377-48-7

MDL:


MFCD00011086

Formula:

Li<sub>2</sub>O<sub>4</sub>S

Formula Weight:

109.94 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form) Powder, Crystals, Crystalline Powder and/or Chunks	Conforms to Requirements
Gravimetric Analysis % Sulfur (With Barium Chloride)	28.4 - 29.9 %
ICP Major Analysis Confirms Lithium and Sulfur Components	Confirmed
Purity > = 99.99% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	$\leq$ 100.0 ppm

Specification: PRD.2.ZQ5.10000022628

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.