

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Human IgG (whole molecule) antibody produced in goat - whole antiserum

Product Number:

I1011

Storage Temperature:

-20 °C

TEST**Specification**

| | |
|--------------------------------|--|
| Appearance (Form) | Liquid |
| Appearance (Turbidity) | Clear to Slightly Hazy |
| Appearance (Color) | Very Light Brown to Very Dark Brown and Dark Brown-Red to Very Dark Brown-Red |
| Protein by Biuret | 40 - 100 mg/ml |
| Antibody Concentration(Method) | Quantitative Precipitin Assay |
| Antibody Concentration | 4.0 - 6.0 mg/ml |
| ELISA - Titer | 1:40,000-1:80,000 |
| Recommended Retest Period | ----- |
| 5 Years | |

Specification: PRE.0.ZQ5.10000020879

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.