

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Hederagenin - ≥97% (HPLC)

Product Number: H3916
CAS Number: 465-99-6
Formula: C₃₀H₄₈O₄
Formula Weight: 472.70 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	White to Light Brown
Appearance (Form)	Powder
Solubility (Color)	Colorless to Light Brown
Solubility (Turbidity)	Clear
1 mg/mL, MeOH	
Water (by Karl Fischer)	≤ 3 %
Purity (HPLC)	≥ 97 %

Specification: PRD.0.ZQ5.10000012162

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.