

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium bromide - $\geq 99.99\%$ trace metals basis

Product Number: 229881
CAS Number: 7647-15-6
MDL: MFCD00003475
Formula: BrNa
Formula Weight: 102.89 g/mol

NaBr

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Titration by AgNO ₃ % Br	$\geq 75.7\%$
ICP Major Analysis Confirms Sodium Component	Confirmed
Trace Metal Analysis	≤ 100.0 ppm
Purity $\geq 99.99\%$ Based On Trace Metals Analysis	Meets Requirements

Specification: PRD.0.ZQ5.10000025850

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.