

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)


## Product Specification

Product Name:

Branched polyethylenimine solution – 2mg/mL aqueous solution, suitable for biomedical research

Product Number:

913375



Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)

Colorless to Light Yellow

Appearance (Form)

Liquid

1H NMR Spectrum

Conforms to Structure

Concentration

1.5 - 2.5 mg/ml

by Residue on Evaporation

Recommended Retest Period

2 Years

Specification: PRD.0.ZQ5.10000103654

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

