

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Barium sulfate – ReagentPlus®, 99%

Product Number: 243353
CAS Number: 7727-43-7
MDL: MFCD00003455
Formula: BaO4S
Formula Weight: 233.39 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
X-Ray Diffraction	Conforms to Structure
Assay (USP)	≥ 98.5 %

Specification: PRD.1.ZQ5.10000012590

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.