

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2,4-DIPROPOXYPHENYLBORONIC ACID, > = 95%

Product Number: 680672
CAS Number: 150145-25-8

Formula: C₁₂H₁₉BO₄
Formula Weight: 238.09 g/mol

TEST	Specification
Appearance (Color)	White to Tan
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder, and/or Chunks	
Proton NMR Spectrum	Conforms to Structure
Titration with NaOH	≥ 95.0 %

Specification: PRD.0.ZQ5.10000039514

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.