

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-p53 antibody produced in sheep - fractionated antiserum, buffered aqueous solution

Product Number: **P4235**

MDL: MFCD00164709

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
Solution in phosphate buffered saline containing 0.08% sodium azide	
Concentration	Conforms
1 mg/mL	
Functional Analysis	Conforms
Western Blot	

Specification: PRD.0.ZQ5.10000053806

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.