


Product Specification

Product Name:
Mevastatin - $\geq 98\%$ (HPLC), powder or crystals

Product Number: **M2537**
CAS Number: 73573-88-3
Formula: C₂₃H₃₄O₅
Formula Weight: 390,51 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Form)- Powder or crystals	Conforms
Appearance (Color)- White to Light Yellow	Conforms
Total Impurities by HPLC	$\leq 2.0 \%$
Any Individual Impurity (HPLC)	$\leq 0.5 \%$
Assay (HPLC) (dried basis)	$\geq 98.0 \%$
Identification: IR-Conforms	Conforms
Identification: UV-Conforms	Conforms
Specific Rotation	265 - 290 °
Loss on Drying	$\leq 0.5 \%$
Residue on Ignition	$\leq 0.1 \%$

Specification: PRE.5.ZQ5.10000014791

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

